HeadaTerm bruksanvisning

Snabb bruksanvisning

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Address: Rauhalantie 70, 07900 Loviisa

Tel: 045 609 9000

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European Authorized Representative

Name: KINGSMEAD SERVICE LIMITED
Add: 145-157 St John Street London EC1V 4PY UK
Contact Person: Edison.V
Tel: 0044-2032399738
Email: office@kingsmead-service.com
Effective Dale: 2018 -05-02 Rev. CE-A/2

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Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Attachment: Annex of Report IEC60601-1-2 Manufacturer's Declaration of the EUT

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT and SYSTEMS

1	Guidance and manufacturer's declaration – electromagnetic emission		
2	The YF-HT-W1 TENS device-HeadTerm is intended for use in the electromagnetic environment specified below. The customer or the user of YF-HT-W1 TENS device-HeadTerm should assure that it is used in such an environment.		
3	Emissions test	Compliance	Electromagnetic environment - guidance
4	RF emissions CISPR 11	Group 1	The YF-HT-W1 TENS device-HeadTerm uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR 11	Class B	The YF-HT-W1 TENS device-HeadTerm is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes.
6	Harmonic emissions IEC 61000-3-2	N/A	
7	Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

Recommended separation distances between portable and mobile RF communications equipment and the YF-HT-W1 TENS device-HeadTerm

The YF-HT-W1 TENS device-HeadTerm is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the YF-HT-W1 TENS device-HeadTerm can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the YF-HT-W1 TENS device-HeadTerm as recommended below, according to the maximum output power of the communications equipment

Rated maximum output of transmitter	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{W} \right] \sqrt{P}$	$d = \left[\frac{3.5}{P} \right] \sqrt{P}$	$d = \left[\frac{2}{P} \right] \sqrt{P}$
0.01	/	0.12	0.23
0.1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity

The YF-HT-W1 TENS device-HeadTerm is intended for use in the electromagnetic environment specified below. The customer or the user of the YF-HT-W1 TENS device-HeadTerm should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U _r (>95 % dip in U _r) for 0.5 cycle 40 % U _r (60 % dip in U _r) for 5 cycles 70 % U _r (30 % dip in U _r) for 25 cycles < 5 % U _r (>95 % dip in U _r) for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the YF-HT-W1 TENS device-HeadTerm requires continued operation during power mains interruptions, it is recommended that the YF-HT-W1 TENS device-HeadTerm be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_r is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity

The YF-HT-W1 TENS device-HeadTerm is intended for use in the electromagnetic environment specified below. The customer or the user of the YF-HT-W1 TENS device-HeadTerm should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the YF-HT-W1 TENS device-HeadTerm including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{W} \right] \sqrt{P}$ $d = \left[\frac{3.5}{P} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{2}{P} \right] \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the YF-HT-W1 TENS device-HeadTerm is used exceeds the applicable RF compliance level above, the YF-HT-W1 TENS device-HeadTerm should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the YF-HT-W1 TENS device-HeadTerm.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.